

**ABRF 2009:**  
**Optimization and Application of Existing and Emerging Biotechnologies**  
February 7-10, 2009 – Memphis, Tennessee

**SATELLITE EDUCATIONAL WORKSHOP PROGRAM**  
**(sw3) Recombinant Protein Laboratory**  
(current as of 12/1/08)

February 6-7, 2009  
8:00 am – 4:00 pm  
St Jude Children's Research Hospital.

**Practical Aspects of Recombinant Protein Expression and Purification**

John Hawes, Miami University (organizer),  
Richard Heath, St. Jude Children's Research Hospital (co-organizer),  
James Bryson, Bristol-Myers Squibb, Preston Hensley, Pfizer,  
Cynthia Kinsland, Cornell University, and Francis Rajamohan, Pfizer

**FRIDAY, FEBRUARY 6, 2009 (DAY 1)**

- |                    |  |
|--------------------|--|
| 7:00 am - 12:00 pm | <b>REGISTRATION OPEN</b>   |
| 7:00 - 8:00 am     | <b>CONTINENTAL BREAKFAST</b> – <i>St Jude, DTRC Room T4</i>  |
| 8:00 - 9:00 am     | <b>INTRODUCTION: Basics of Protein Expression</b> – <i>St. Jude, E1004</i><br><b>A. Choice of a Host System</b> <ul style="list-style-type: none"><li>1. Bacteria, yeast, insect and mammalian culture</li><li>2. Equipment needed for core lab</li></ul> <b>B. Biology of Heterologous Protein Overexpression in Bacteria</b> <ul style="list-style-type: none"><li>1. Overview of the pET system</li><li>2. The biggest failure in bacteria: insoluble protein</li><li>3. Tips for shifting the balance to soluble expression</li></ul> <b>C. Details of Lab Work</b> <ul style="list-style-type: none"><li>1. Growth of <i>E. coli</i> in shake flasks</li><li>2. Introduction to 10 L fermentors</li></ul> |
| 9:15 am - 12:00 pm | <b>LAB WORK: Grow and Induce Recombinant Cultures</b> – <i>St. Jude, Lab E8066</i> <ul style="list-style-type: none"><li>1. Hands on lab work in groups of 3</li><li>2. Demonstration of 10 L bench top fermentor</li></ul>  |
| 12:00 - 1:00 pm    | <b>LUNCH</b> – <i>St. Jude, E1004</i>  |
| 1:00 – 1:30 pm     | <b>SOLUBILITY: Screening For Suitable Buffers</b> – <i>St. Jude, E1004</i><br><b>A. Use of Prefilled 96 Well Trays to Select Optimum Buffers for Purified Proteins</b><br><b>B. Details of Lab Work</b>  |
| 1:30 – 4:00 pm     | <b>LAB WORK</b> – <i>St. Jude, Lab E8066</i><br><b>A. Solubility Screening</b> <ul style="list-style-type: none"><li>1. Set up solubility screens using supplied purified protein</li></ul> <b>B. Growth and Harvest of Cultures</b> <ul style="list-style-type: none"><li>1. Continuation of morning growths</li><li>2. Harvest by centrifugation</li></ul>   |

**SATURDAY, FEBRUARY 7, 2009 (DAY 2)**

7:00 am - 12:00 pm      **REGISTRATION OPEN**

7:00 - 8:00 am      **CONTINENTAL BREAKFAST** – *St. Jude, E1004*

8:00 - 9:00 am      **BASICS AND TIPS OF AFFINITY PURIFICATION**  
1. Introduction to metal chelation affinity chromatography  
2. Other chromatography methods  
3. Benchtop versus AKTA

9:15 am - 12:00 pm      **LAB WORK** – *St. Jude, Lab E8066*  
**A. Cell Lysis**  
1. Mechanical (demonstration)  
2. Detergent-based  
  
**B. Protein Purification**  
1. Running a HisTrap column on an AKTA (demonstration)  
2. Bench-top purification using MCAC resins

12:00 - 1:00 pm      **LUNCH** – *St. Jude, E1004*

1:00 - 1:30 pm      **ASK THE EXPERTS SESSION** – *St. Jude, E1004*  
Round-table discussion on protein expression/purification as a shared resource

1:30 - 4:00 pm      **LAB WORK** – *St. Jude, Lab E8066*  
**A. SDS-PAGE Analysis of Purified Protein**  
**B. Biochemical Assay of Activity**